

产品名称: KCTD3 Rabbit Monoclonal Antibody
产品货号: AMRe02188



产品概述 (Summary)

产品名称 (Production Name)	KCTD3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KCTD3
别名 (Alternative Names)	Renal carcinoma antigen NY-REN-45
基因 ID (Gene ID)	51133
蛋白 ID (SwissProt ID)	Q9Y597.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 89 kDa; Observed MW: 89 kDa

研究背景 (Background)

Accessory subunit of potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 3 (HCN3) up-

产品名称: KCTD3 Rabbit Monoclonal Antibody
产品货号: AMRe02188




regulating its cell-surface expression and current density without affecting its voltage dependence and kinetics.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of KCTD3 in K562 lysates using KCTD3 antibody.

注意事项 (Note)

For research use only .