

产品名称: KCTD9 Rabbit Monoclonal Antibody
产品货号: AMRe02189



产品概述 (Summary)

产品名称 (Production Name)	KCTD9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KCTD9
别名 (Alternative Names)	BTB/POZ domain-containing protein KCTD9
基因 ID (Gene ID)	54793
蛋白 ID (SwissProt ID)	Q7L273.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43 kDa

研究背景 (Background)

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex, which mediates the

产品名称: KCTD9 Rabbit Monoclonal Antibody
产品货号: AMRe02189




ubiquitination of target proteins, leading to their degradation by the proteasome.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of KCTD9 in K562, HeLa lysates using KCTD9 antibody.



Immunocytochemistry analysis of KCTD9(green) in HeLa using KCTD9 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .