

产品名称: KIFC1 Rabbit Monoclonal Antibody
产品货号: AMRe02194



产品概述 (Summary)

产品名称 (Production Name)	KIFC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KIFC1
别名 (Alternative Names)	HSET; KNSL2
基因 ID (Gene ID)	3833
蛋白 ID (SwissProt ID)	Q9BW19.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 74 kDa; Observed MW: 74 kDa

研究背景 (Background)

Minus end-directed microtubule-dependent motor required for bipolar spindle formation. May contribute to movement of

产品名称: KIFC1 Rabbit Monoclonal Antibody
产品货号: AMRe02194




early endocytic vesicles.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of KIFC1 in K562, 3T3, HeLa lysates using KIFC1 antibody.



Immunocytochemistry analysis of KIFC1(green) in HeLa using KIFC1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using KIFC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: KIFC1 Rabbit Monoclonal Antibody
产品货号: AMRe02194



注意事项 (Note)

For research use only .