

产品名称: L1CAM Rabbit Monoclonal Antibody
产品货号: AMRe02198



产品概述 (Summary)

产品名称 (Production Name)	L1CAM Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	L1CAM
别名 (Alternative Names)	L1CAM; CAML1; MIC5; Neural cell adhesion molecule L1; N-CAM-L1; NCAM-L1; CD antigen CD171
基因 ID (Gene ID)	3897
蛋白 ID (SwissProt ID)	P32004.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 140 kDa; Observed MW: 140 kDa

研究背景 (Background)

产品名称: L1CAM Rabbit Monoclonal Antibody
产品货号: AMRe02198




Cell adhesion molecule with an important role in the development of the nervous system. Involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc. Binds to axonin on neurons.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of L1CAM in A375 lysates using L1CAM antibody.

注意事项 (Note)

For research use only .