

产品名称: Lamin B2 Rabbit Monoclonal Antibody
产品货号: AMRe02203



产品概述 (Summary)

产品名称 (Production Name)	Lamin B2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LMNB2
别名 (Alternative Names)	LMN2; LAMB2
基因 ID (Gene ID)	84823
蛋白 ID (SwissProt ID)	Q03252.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

Lamins are components of the nuclear lamina, a fibrous layer on the nucleoplasmic side of the inner nuclear membrane,

产品名称: Lamin B2 Rabbit Monoclonal Antibody
产品货号: AMRe02203




which is thought to provide a framework for the nuclear envelope and may also interact with chromatin.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of Lamin B2 in C6, 3T3, HeLa lysates using Lamin B2 antibody.



Immunocytochemistry analysis of Lamin B2(green) in HeLa using Lamin B2 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Lamin B2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Lamin B2 Rabbit Monoclonal Antibody
产品货号: AMRe02203



注意事项 (Note)

For research use only .