

产品名称: LGALS14 Rabbit Monoclonal Antibody
产品货号: AMRe02212



产品概述 (Summary)

产品名称 (Production Name)	LGALS14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LGALS14
别名 (Alternative Names)	Charcot-Leyden crystal protein 2; CLC2; Galectin-14; Gal-14
基因 ID (Gene ID)	56891
蛋白 ID (SwissProt ID)	Q8TCE9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 16 kDa

研究背景 (Background)

Binds beta-galactoside and lactose.


产品名称: LGALS14 Rabbit Monoclonal Antibody
产品货号: AMRe02212




研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Immunocytochemistry analysis of LGALS14 (green) in K562 using LGALS14 antibody, and DAPI (blue).



Western blot analysis of LGALS14 in C6 lysates using LGALS14 antibody.

注意事项 (Note)

For research use only .