

产品名称: LRRC15 Rabbit Monoclonal Antibody
产品货号: AMRe02221



产品概述 (Summary)

产品名称 (Production Name)	LRRC15 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LRRC15
别名 (Alternative Names)	LIB
基因 ID (Gene ID)	131578
蛋白 ID (SwissProt ID)	Q8TF66.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 64 kDa; Observed MW: 64 kDa

研究背景 (Background)

Enables collagen binding activity; fibronectin binding activity; and laminin binding activity. Involved in negative regulation

产品名称: LRRC15 Rabbit Monoclonal Antibody
产品货号: AMRe02221




of protein localization to plasma membrane; positive regulation of cell migration; and receptor-mediated virion attachment to host cell. Located in extracellular exosome. Colocalizes with collagen-containing extracellular matrix.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of LRRC15 in rat Brain, C6 lysates using LRRC15 antibody.

注意事项 (Note)

For research use only .