

产品名称: Lumican Rabbit Monoclonal Antibody  
产品货号: AMRe02222



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Lumican Rabbit Monoclonal Antibody     |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF                          |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                           |
|------------------------|---------------------------|
| 基因名 (Gene Name)        | LUM                       |
| 别名 (Alternative Names) | LDC; LUM; Lumican; SLRR2D |
| 基因 ID (Gene ID)        | 4060                      |
| 蛋白 ID (SwissProt ID)   | P51884.                   |

## 产品应用 (Application)

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 38 kDa; Observed MW: 60 kDa       |

## 研究背景 (Background)

Belongs to the small leucine-rich proteoglycan (SLRP) family. SLRP class II subfamily. Contains 11 LRR (leucine-rich) repeats.

产品名称: Lumican Rabbit Monoclonal Antibody  
产品货号: AMRe02222




Contains 1 LRRNT domain.


## 研究领域 (Research Area)

Cancer


## 图片 (Image Data)



Western blot analysis of Lumican in K562, C6, 3T3, HeLa lysates using Lumican antibody.



Immunocytochemistry analysis of Lumican(green) in HeLa using Lumican antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Lumican antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

**产品名称: Lumican Rabbit Monoclonal Antibody**  
**产品货号: AMRe02222**

---



**注意事项 (Note)**

For research use only .