

产品名称: LYVE1 Rabbit Monoclonal Antibody
产品货号: AMRe02224



产品概述 (Summary)

产品名称 (Production Name)	LYVE1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lyve1
别名 (Alternative Names)	CRSBP-1; Extracellular link domain-containing protein 1
基因 ID (Gene ID)	114332.0
蛋白 ID (SwissProt ID)	Q8BHC0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa; Observed MW: 35-70 kDa

研究背景 (Background)

Ligand-specific transporter trafficking between intracellular organelles (TGN) and the plasma membrane. Plays a role in

产品名称: LYVE1 Rabbit Monoclonal Antibody
产品货号: AMRe02224




autocrine regulation of cell growth mediated by growth regulators containing cell surface retention sequence binding (CRS).


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse lung using LYVE1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of LYVE1 in mouse Lung lysates using LYVE1 antibody.

注意事项 (Note)

For research use only .