

产品名称: MAGEA3 Rabbit Monoclonal Antibody
产品货号: AMRe02226



产品概述 (Summary)

产品名称 (Production Name)	MAGEA3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAGEA3
别名 (Alternative Names)	HIP8; HYPD; CT1.3; MAGE3; MAGEA6
基因 ID (Gene ID)	4102
蛋白 ID (SwissProt ID)	P43357.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa; Observed MW: 45 kDa

研究背景 (Background)

Proposed to enhance ubiquitin ligase activity of RING-type zinc finger-containing E3 ubiquitin-protein ligases. May

产品名称: MAGEA3 Rabbit Monoclonal Antibody
产品货号: AMRe02226




enhance ubiquitin ligase activity of TRIM28 and stimulate p53/TP53 ubiquitination by TRIM28. Proposed to act through recruitment and/or stabilization of the Ubl-conjugating enzyme (E2) at the E3:substrate complex. May play a role in embryonal development and tumor transformation or aspects of tumor progression. In vitro promotes cell viability in melanoma cell lines. Antigen recognized on a melanoma by autologous cytolytic T-lymphocytes.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of MAGEA3 in HeLa lysates using MAGEA3 antibody.

注意事项 (Note)

For research use only .