

产品名称: c-Met Rabbit Monoclonal Antibody
产品货号: AMRe02248



产品概述 (Summary)

产品名称 (Production Name)	c-Met Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MET MET; Hepatocyte growth factor receptor; HGF receptor; HGF/SF receptor;
别名 (Alternative Names)	Proto-oncogene c-Met; Scatter factor receptor; SF receptor; Tyrosine-protein kinase Met
基因 ID (Gene ID)	4233
蛋白 ID (SwissProt ID)	P08581.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 156 kDa; Observed MW: 170 kDa

产品名称: c-Met Rabbit Monoclonal Antibody
产品货号: AMRe02248




研究背景 (Background)

The proto-oncogene MET product is the hepatocyte growth factor receptor and encodes tyrosine-kinase activity. The primary single chain precursor protein is post-translationally cleaved to produce the alpha and beta subunits, which are disulfide linked to form the mature receptor.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Met (c-Met) in K562, C6, 3T3 lysates using c-Met antibody.



Immunocytochemistry analysis of Met (c-Met)(green) in HeLa using Met (c-Met) antibody, and DAPI(blue)



产品名称: c-Met Rabbit Monoclonal Antibody
产品货号: AMRe02248



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Met (c-Met) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .