

产品名称: MIB1 Rabbit Monoclonal Antibody
产品货号: AMRe02254



产品概述 (Summary)

产品名称 (Production Name)	MIB1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MIB1
别名 (Alternative Names)	Dip1; LVNC7; MIB; mib1; ZZANK2; ZZZ6; DAPK-interacting protein 1; Mind bomb homolog 1; RING-type E3 ubiquitin transferase MIB1
基因 ID (Gene ID)	57534
蛋白 ID (SwissProt ID)	Q86YT6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 110 kDa; Observed MW: 110 kDa

研究背景 (Background)

产品名称: MIB1 Rabbit Monoclonal Antibody
产品货号: AMRe02254




E3 ubiquitin-protein ligase that mediates ubiquitination of Delta receptors, which act as ligands of Notch proteins. Positively regulates the Delta-mediated Notch signaling by ubiquitinating the intracellular domain of Delta, leading to endocytosis of Delta receptors.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Mib1/Mindbomb in K562, Hela lysates using MIB1 antibody.

注意事项 (Note)

For research use only .