

产品名称: MMP12 Rabbit Monoclonal Antibody
产品货号: AMRe02263



产品概述 (Summary)

产品名称 (Production Name)	MMP12 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MMP12
别名 (Alternative Names)	Macrophage metalloelastase; MME; 3.4.24.65; Macrophage elastase; ME; hME; Matrix metalloproteinase-12; MMP-12; MMP12; HME
基因 ID (Gene ID)	4321
蛋白 ID (SwissProt ID)	P39900.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54,45 kDa

研究背景 (Background)

产品名称: MMP12 Rabbit Monoclonal Antibody
产品货号: AMRe02263




MMP-12 lacks a transmembrane domain and furin cleavage site. The zymogen for MMP-12 is about 54 kDa, and is quickly activated to the 45 kDa form, and this breaks down to cascade of active forms, ending with the common 22 kDa form. Stimulated macrophages produce MMP-12; it has also been found in osteosarcoma cells, synovial fibroblasts and lung fibroblasts.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of MMP12 in HeLa lysates using MMP12 antibody.



Western blot analysis of MMP12 in CHO-K1, C6 lysates using MMP12 antibody.

注意事项 (Note)

For research use only .