

产品名称: Myelin Basic Protein Rabbit Monoclonal Antibody
产品货号: AMRe02294



产品概述 (Summary)

产品名称 (Production Name)	Myelin Basic Protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MBP
别名 (Alternative Names)	Myelin membrane encephalitogenic protein; Myelin A1 protein
基因 ID (Gene ID)	4155
蛋白 ID (SwissProt ID)	P02686.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 18 kDa

研究背景 (Background)

MBP belongs to the myelin basic protein family. The classic group of MBP isoforms (isoform 4-isoform 14) are the most

产品名称: Myelin Basic Protein Rabbit Monoclonal Antibody
产品货号: AMRe02294




abundant protein components of the myelin membrane in the CNS. They have a role in both its formation and stabilization. The smaller isoforms might have an important role in remyelination of denuded axons in multiple sclerosis.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Myelin Basic Protein in HeLa, CHO-K1, C6 lysates using Myelin Basic Protein antibody.

注意事项 (Note)

For research use only .