

产品名称: NCK1 Rabbit Monoclonal Antibody
产品货号: AMRe02309



产品概述 (Summary)

产品名称 (Production Name)	NCK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCK1
别名 (Alternative Names)	NCK; nck-1; NCKalpha
基因 ID (Gene ID)	4690
蛋白 ID (SwissProt ID)	P16333.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43 kDa

研究背景 (Background)

Adapter protein which associates with tyrosine-phosphorylated growth factor receptors or their cellular substrates.

产品名称: NCK1 Rabbit Monoclonal Antibody
产品货号: AMRe02309




Maintains low levels of EIF2S1 phosphorylation by promoting its dephosphorylation by PP1.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of NCK1 (green) in HeLa using NCK1 antibody, and DAPI(blue).



Western blot analysis of Nck in K562, rat Brain, C6, 3T3, HeLa lysates using Nck antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Nck antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: NCK1 Rabbit Monoclonal Antibody
产品货号: AMRe02309



注意事项 (Note)

For research use only .