

产品名称: NDRG2 Rabbit Monoclonal Antibody
产品货号: AMRe02310



产品概述 (Summary)

产品名称 (Production Name)	NDRG2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDRG2
别名 (Alternative Names)	SYLD
基因 ID (Gene ID)	57447
蛋白 ID (SwissProt ID)	Q9UN36.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)

Contributes to the regulation of the Wnt signaling pathway. Down-regulates CTNNB1-mediated transcriptional activation

产品名称: NDRG2 Rabbit Monoclonal Antibody
产品货号: AMRe02310




of target genes, such as CCND1, and may thereby act as tumor suppressor. May be involved in dendritic cell and neuron differentiation.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of NDRG2 in rat brain lysates using NDRG2 antibody.



Western blot analysis of NDRG2 in mouse brain lysates using NDRG2 antibody.

注意事项 (Note)

For research use only .