

产品名称: NDUFAF4 Rabbit Monoclonal Antibody
产品货号: AMRe02314



产品概述 (Summary)

产品名称 (Production Name)	NDUFAF4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFAF4
别名 (Alternative Names)	My013; C6orf66; HRPAP20; HSPC125; bA22L21.1
基因 ID (Gene ID)	29078
蛋白 ID (SwissProt ID)	Q9P032.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 20 kDa; Observed MW: 20 kDa

研究背景 (Background)

Involved in the assembly of mitochondrial NADH:ubiquinone oxidoreductase complex (complex I). May be involved in cell

产品名称: NDUFAF4 Rabbit Monoclonal Antibody
产品货号: AMRe02314




proliferation and survival of hormone-dependent tumor cells. May be a regulator of breast tumor cell invasion.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Immunocytochemistry analysis of NDUFAF4 (green) in HCT116 using NDUFAF4 antibody, and DAPI(blue).



Western blot analysis of NDUFAF4 in K562, HeLa lysates using NDUFAF4 antibody.

注意事项 (Note)

For research use only .