

产品名称: NDUFB4 Rabbit Monoclonal Antibody
产品货号: AMRe02319



产品概述 (Summary)

产品名称 (Production Name)	NDUFB4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFB4
别名 (Alternative Names)	B15; CI-B15
基因 ID (Gene ID)	4710
蛋白 ID (SwissProt ID)	O95168.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

研究背景 (Background)

Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not

产品名称: NDUFB4 Rabbit Monoclonal Antibody
产品货号: AMRe02319




to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.


研究领域 (Research Area)

Tags & Cell Markers


图片 (Image Data)



Immunocytochemistry analysis of NDUFB4 (green) in Hela using NDUFB4 antibody, and DAPI (blue).




Western blot analysis of NDUFB4 in mouse skeletal muscle lysates using NDUFB4 antibody.



Western blot analysis of NDUFB4 in K562, rat Brain lysates using NDUFB4 antibody

产品名称: NDUFB4 Rabbit Monoclonal Antibody
产品货号: AMRe02319



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using NDUFB4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .