

产品名称: NDUFS6 Rabbit Monoclonal Antibody
产品货号: AMRe02324



产品概述 (Summary)

产品名称 (Production Name)	NDUFS6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFS6
别名 (Alternative Names)	CI-13kA; CI13KDA; CI-13kD-A
基因 ID (Gene ID)	4726
蛋白 ID (SwissProt ID)	O75380.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 14 kDa; Observed MW: 14 kDa

研究背景 (Background)

Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not

产品名称: NDUFS6 Rabbit Monoclonal Antibody
产品货号: AMRe02324




to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.


研究领域 (Research Area)

Tags & Cell Markers


图片 (Image Data)



Western blot analysis of NDUFS6 in mouse heart lysates using NDUFS6 antibody.



Western blot analysis of NDUFS6 in HeLa lysates using NDUFS6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using NDUFS6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: NDUFS6 Rabbit Monoclonal Antibody
产品货号: AMRe02324



注意事项 (Note)

For research use only .