

产品名称: NEAS Rabbit Monoclonal Antibody
产品货号: AMRe02326



产品概述 (Summary)

产品名称 (Production Name)	NEAS Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SPTAN1
别名 (Alternative Names)	SPTAN1; NEAS; SPTA2; Spectrin alpha chain; non-erythrocytic 1; Alpha-II spectrin; Fodrin alpha chain; Spectrin; non-erythroid alpha subunit
基因 ID (Gene ID)	6709
蛋白 ID (SwissProt ID)	Q13813.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 285 kDa; Observed MW: 285 kDa

研究背景 (Background)

产品名称: NEAS Rabbit Monoclonal Antibody
产品货号: AMRe02326




Fodrin (also named nonerythroid spectrin) is a universally expressed membrane-associated cytoskeletal protein consisting of alpha- and beta-subunits (1). This protein is important for maintaining normal membrane structure and supporting cell surface protein function (1). Alpha-fodrin is one of the primary targets cleaved by caspases during apoptosis.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of NEAS in Hela, Hacat lysates using NEAS antibody.



Western blot analysis of NEAS in mouse brain, rat brain lysates using NEAS antibody.

注意事项 (Note)

For research use only .