

产品名称: NEK4 Rabbit Monoclonal Antibody
产品货号: AMRe02327



产品概述 (Summary)

产品名称 (Production Name)	NEK4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEK4
别名 (Alternative Names)	NRK2; STK2; pp12301
基因 ID (Gene ID)	6787
蛋白 ID (SwissProt ID)	P51957.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 95 kDa; Observed MW: 95 kDa

研究背景 (Background)

Protein kinase that seems to act exclusively upon threonine residues . Required for normal entry into proliferative arrest

产品名称: NEK4 Rabbit Monoclonal Antibody
产品货号: AMRe02327




after a limited number of cell divisions, also called replicative senescence. Required for normal cell cycle arrest in response to double-stranded DNA damage.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of NEK4 in Jurkat, HepG2, HeLa lysates using NEK4 antibody.

注意事项 (Note)

For research use only .