

产品名称: NeuN Rabbit Monoclonal Antibody
产品货号: AMRe02334



产品概述 (Summary)

产品名称 (Production Name)	NeuN Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBFOX3
别名 (Alternative Names)	FLJ56884; FLJ58356; FOX3; HRNBP3; RBFOX3; NEUN; FOX-3
基因 ID (Gene ID)	146713
蛋白 ID (SwissProt ID)	A6NFN3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 46-55 kDa

研究背景 (Background)

NeuN, also named as FOX3 and RBFOX3, is neuronal-specific nuclear protein with MW 38-50kd(~48kd). It is one of a family

产品名称: NeuN Rabbit Monoclonal Antibody
产品货号: AMRe02334




of 3 mammalian Fox homologues. FOX3(Or NeuN) functions in RNA-binding protein that regulates alternative splicing events.


研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Immunocytochemistry analysis of NeuN (green) in SH-SY5Y using NeuN antibody, and DAPI (blue).



Western blot analysis of NeuN in mouse brain lysates using NeuN antibody.



Western blot analysis of NeuN in rat Brain lysates using NeuN antibody

产品名称: NeuN Rabbit Monoclonal Antibody
产品货号: AMRe02334



注意事项 (Note)

For research use only .