

产品名称: Niemann Pick C2 Rabbit Monoclonal Antibody
产品货号: AMRe02344



产品概述 (Summary)

产品名称 (Production Name)	Niemann Pick C2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPC2
别名 (Alternative Names)	HE1; EDDM1
基因 ID (Gene ID)	10577
蛋白 ID (SwissProt ID)	P61916.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

研究背景 (Background)

This gene encodes a protein containing a lipid recognition domain. The encoded protein may function in regulating the

产品名称: Niemann Pick C2 Rabbit Monoclonal Antibody
产品货号: AMRe02344




transport of cholesterol through the late endosomal/lysosomal system. Mutations in this gene have been associated with Niemann-Pick disease, type C2 and frontal lobe atrophy. [provided by RefSeq, Jul 2008]


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Niemann Pick C2 in C6, 3T3, HeLa lysates using Niemann Pick C2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Niemann Pick C2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .