

产品名称: THO Complex Subunit 1 Rabbit Monoclonal Antibody
产品货号: AMRe02359



产品概述 (Summary)

产品名称 (Production Name)	THO Complex Subunit 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THOC1
别名 (Alternative Names)	Death domain containing protein p84N5; HPR1; hTRES84; Nuclear matrix protein p84; p84N5; THO complex 1; Tho1; Thoc1
基因 ID (Gene ID)	9984
蛋白 ID (SwissProt ID)	Q96FV9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa; Observed MW: 84 kDa

产品名称: THO Complex Subunit 1 Rabbit Monoclonal Antibody
产品货号: AMRe02359




研究背景 (Background)

Regulates transcriptional elongation of a subset of genes. Participates in an apoptotic pathway which is characterized by activation of caspase-6, increases in the expression of BAK1 and BCL2L1 and activation of NF-kappa-B.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Nuclear Matrix Protein p84 in Jurkat, HeLa, COS7 lysates using THO Complex Subunit 1 antibody.



Immunocytochemistry analysis of Nuclear Matrix Protein p84(green) in HeLa using Nuclear Matrix Protein p84 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .