

产品名称: Nuf2 Rabbit Monoclonal Antibody
产品货号: AMRe02366



产品概述 (Summary)

产品名称 (Production Name)	Nuf2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUF2
别名 (Alternative Names)	CDCA1; CT106; hNuf2; hNuf2R; hsNuf2; NUF2; NUF2R
基因 ID (Gene ID)	83540
蛋白 ID (SwissProt ID)	Q9BZD4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54 kDa

研究背景 (Background)

Acts as a component of the essential kinetochore-associated NDC80 complex, which is required for chromosome

产品名称: Nuf2 Rabbit Monoclonal Antibody
产品货号: AMRe02366




segregation and spindle checkpoint activity. Required for kinetochore integrity and the organization of stable microtubule binding sites in the outer plate of the kinetochore.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Nuf2 in K562 lysates using Nuf2 antibody.

注意事项 (Note)

For research use only .