

产品名称: NXF1 Rabbit Monoclonal Antibody
产品货号: AMRe02370



产品概述 (Summary)

产品名称 (Production Name)	NXF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NXF1
别名 (Alternative Names)	DmNXF1; MEX67; Mex67p; Mvb1; NXF1; Sbr; TAP
基因 ID (Gene ID)	10482
蛋白 ID (SwissProt ID)	Q9UBU9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

Involved in the nuclear export of mRNA species bearing retroviral constitutive transport elements (CTE) and in the export of

产品名称: NXF1 Rabbit Monoclonal Antibody
产品货号: AMRe02370




mRNA from the nucleus to the cytoplasm. The NXF1-NXT1 heterodimer is involved in the export of HSP70 mRNA in conjunction with THOC4 and THOC5.


研究领域 (Research Area)

Tags & Cell Markers


图片 (Image Data)



Immunocytochemistry analysis of NXF1 (green) in Jurkat using NXF1 antibody, and DAPI (blue).



Western blot analysis of NXF1 in K562 lysates using NXF1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using NXF1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: NXF1 Rabbit Monoclonal Antibody
产品货号: AMRe02370



注意事项 (Note)

For research use only .