

产品名称: Oct-1 Rabbit Monoclonal Antibody
产品货号: AMRe02371



产品概述 (Summary)

产品名称 (Production Name)	Oct-1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POU2F1 POU2F1; OCT1; OTF1; POU domain; class 2; transcription factor 1; NF-A1;
别名 (Alternative Names)	Octamer-binding protein 1; Oct-1; Octamer-binding transcription factor 1; OTF-1; POU2F2; OCT2; OTF2; POU domain; class 2; transcription factor 2; Lymphoid-restricted im
基因 ID (Gene ID)	5451
蛋白 ID (SwissProt ID)	P14859.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa; Observed MW: 90 kDa

产品名称: Oct-1 Rabbit Monoclonal Antibody
产品货号: AMRe02371




研究背景 (Background)

Transcription factor that binds to the octamer motif (5'-ATTGTCAT-3') and activates the promoters of the genes for some small nuclear RNAs (snRNA) and of genes such as those for histone H2B and immunoglobulins.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Oct1 in HeLa, 293 lysates using Oct1 antibody.

注意事项 (Note)

For research use only .