

产品名称: ORP1 Rabbit Monoclonal Antibody
产品货号: AMRe02374



产品概述 (Summary)

产品名称 (Production Name)	ORP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OSBPL1A
别名 (Alternative Names)	ORP1; OSBP8; OSBPL1; OSBPL1B
基因 ID (Gene ID)	114876
蛋白 ID (SwissProt ID)	Q9BXW6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 108 kDa; Observed MW: 109,50 kDa

研究背景 (Background)

Binds phospholipids; exhibits strong binding to phosphatidic acid and weak binding to phosphatidylinositol 3-phosphate

产品名称: ORP1 Rabbit Monoclonal Antibody
产品货号: AMRe02374




(By similarity).


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of ORP1 (green) in 293t using ORP1 antibody, and DAPI(blue).



Western blot analysis of ORP1 in mouse brain, rat testis lysates using ORP1 antibody.



Western blot analysis of ORP1 in rat Brain, C6, Hela lysates using ORP1 antibody

产品名称: ORP1 Rabbit Monoclonal Antibody
产品货号: AMRe02374



注意事项 (Note)

For research use only .