

产品名称: ORP150 Rabbit Monoclonal Antibody
产品货号: AMRe02376



产品概述 (Summary)

产品名称 (Production Name)	ORP150 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HYOU1
别名 (Alternative Names)	Grp170; HSP12A; Hyou1; Orp150
基因 ID (Gene ID)	10525
蛋白 ID (SwissProt ID)	Q9Y4L1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 111 kDa; Observed MW: 150 kDa

研究背景 (Background)

Has a pivotal role in cytoprotective cellular mechanisms triggered by oxygen deprivation. May play a role as a molecular

产品名称: ORP150 Rabbit Monoclonal Antibody
产品货号: AMRe02376




chaperone and participate in protein folding.


研究领域 (Research Area)

Cardiovascular


图片 (Image Data)



Western blot analysis of ORP150 in HeLa lysates using ORP150 antibody.



Immunocytochemistry analysis of ORP150(green) in HeLa using ORP150 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using ORP150 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: ORP150 Rabbit Monoclonal Antibody
产品货号: AMRe02376



注意事项 (Note)

For research use only .