

产品名称: p73 Rabbit Monoclonal Antibody
产品货号: AMRe02390



产品概述 (Summary)

产品名称 (Production Name)	p73 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP,ChIP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TP73
别名 (Alternative Names)	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein
基因 ID (Gene ID)	7161
蛋白 ID (SwissProt ID)	O15350.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50,ChIP 1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 80,70 kDa

研究背景 (Background)

产品名称: p73 Rabbit Monoclonal Antibody
产品货号: AMRe02390




p73 IS a protein of the p53 family of proteins. Participates in the apoptotic response to DNA damage. Phosphorylated in a cell cycle-dependent manner and negatively regulated by CDKs. When overproduced, activates transcription from p53-responsive promoters and induces apoptosis. May be a tumor suppressor protein. Seven alternatively-spliced isoforms have been described.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of p73 in 293 lysates using p73 antibody.

注意事项 (Note)

For research use only .