

产品名称: PRKRA Rabbit Monoclonal Antibody
产品货号: AMRe02393



产品概述 (Summary)

产品名称 (Production Name)	PRKRA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRKRA
别名 (Alternative Names)	RAX; PACT; DYT16; HSD14
基因 ID (Gene ID)	8575
蛋白 ID (SwissProt ID)	O75569.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 34 kDa

研究背景 (Background)

Activates EIF2AK2/PKR in the absence of double stranded RNA (dsRNA), leading to phosphorylation of EIF2S1/EIF2- α

产品名称: PRKRA Rabbit Monoclonal Antibody
产品货号: AMRe02393




and inhibition of translation and induction of apoptosis.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of PACT (Phospho- KR activating protein)/PRKRA in C6, 3T3, HeLa lysates using PRKRA antibody.

注意事项 (Note)

For research use only .