

产品名称: PCGF6 Rabbit Monoclonal Antibody
产品货号: AMRe02406



产品概述 (Summary)

产品名称 (Production Name)	PCGF6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Pcgf6
别名 (Alternative Names)	RING finger protein 134; Mel18 and Bmi1-like RING finger
基因 ID (Gene ID)	71041.0
蛋白 ID (SwissProt ID)	Q99NA9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 40 kDa

研究背景 (Background)

Transcriptional repressor. May modulate the levels of histone H3K4Me3 by activating KDM5D histone demethylase.

产品名称: PCGF6 Rabbit Monoclonal Antibody
产品货号: AMRe02406




Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of PCGF6 in PC-12 lysates using PCGF6 antibody.

注意事项 (Note)

For research use only .