

产品名称: PCSK9 Rabbit Monoclonal Antibody
产品货号: AMRe02409



产品概述 (Summary)

产品名称 (Production Name)	PCSK9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PCSK9
别名 (Alternative Names)	HCHOLA3; LDLCQ1; NARC1; Neural apoptosis regulated convertase 1; PC9; PCSK9; Proprotein convertase 9
基因 ID (Gene ID)	255738
蛋白 ID (SwissProt ID)	Q8NBP7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 74 kDa; Observed MW: 65,80 kDa

研究背景 (Background)

产品名称: PCSK9 Rabbit Monoclonal Antibody
产品货号: AMRe02409




May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of PCSK9 in HeLa lysates using PCSK9 antibody.

注意事项 (Note)

For research use only .