

产品名称: PEX11B Rabbit Monoclonal Antibody
产品货号: AMRe02427



产品概述 (Summary)

产品名称 (Production Name)	PEX11B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PEX11B
别名 (Alternative Names)	PEX14B; PEX11-BETA
基因 ID (Gene ID)	8799
蛋白 ID (SwissProt ID)	O96011.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

研究背景 (Background)

Involved in peroxisomal proliferation (PubMed:9792670). May regulate peroxisome division by recruiting the dynamin-

产品名称: PEX11B Rabbit Monoclonal Antibody
产品货号: AMRe02427




related GTPase DNM1L to the peroxisomal membrane (PubMed:12618434). Promotes membrane protrusion and elongation on the peroxisomal surface (PubMed:20826455).


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of PEX11B in T47D, HL-60 lysates using PEX11B antibody.



Western blot analysis of PEX11B in 3T3 lysates using PEX11B antibody.

注意事项 (Note)

For research use only .