

产品名称: PGM1 Rabbit Monoclonal Antibody
产品货号: AMRe02431



产品概述 (Summary)

产品名称 (Production Name)	PGM1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PGM1
别名 (Alternative Names)	Phosphoglucomutase-1; PGM 1; Glucose phosphomutase 1
基因 ID (Gene ID)	5236
蛋白 ID (SwissProt ID)	P36871.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 61 kDa; Observed MW: 61 kDa

研究背景 (Background)

This enzyme participates in both the breakdown and synthesis of glucose.


产品名称: PGM1 Rabbit Monoclonal Antibody
产品货号: AMRe02431




研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of PGM1 in C6, 3T3, HeLa lysates using PGM1 antibody.



Immunocytochemistry analysis of PGM1(green) in HeLa using PGM1 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .