

产品名称: PICK1 Rabbit Monoclonal Antibody
产品货号: AMRe02440



产品概述 (Summary)

产品名称 (Production Name)	PICK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PICK1
别名 (Alternative Names)	PICK 1; Pick1; PRKCABP
基因 ID (Gene ID)	9463
蛋白 ID (SwissProt ID)	Q9NRD5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 50 kDa

研究背景 (Background)

Probable adapter protein that bind to and organize the subcellular localization of a variety of membrane proteins

产品名称: PICK1 Rabbit Monoclonal Antibody
产品货号: AMRe02440




containing some PDZ recognition sequence. Involved in the clustering of various receptors, possibly by acting at the receptor internalization level. Plays a role in synaptic plasticity by regulating the trafficking and internalization of AMPA receptors.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of PICK1 in K562, HeLa lysates using PICK1 antibody.

注意事项 (Note)

For research use only .