

产品名称: PMCA1 Rabbit Monoclonal Antibody
产品货号: AMRe02454



产品概述 (Summary)

产品名称 (Production Name)	PMCA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATP2B1
别名 (Alternative Names)	ATP2B1; PMCA1
基因 ID (Gene ID)	490
蛋白 ID (SwissProt ID)	P20020.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 135 kDa; Observed MW: 130-250 kDa

研究背景 (Background)

This magnesium-dependent enzyme catalyzes the hydrolysis of ATP coupled with the transport of calcium out of the cell.


产品名称: PMCA1 Rabbit Monoclonal Antibody
产品货号: AMRe02454




研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of PMCA1 in mouse brain, rat lung lysates using PMCA1 antibody.



Western blot analysis of PMCA1 in U87-MG, Human placenta lysates using PMCA1 antibody.

注意事项 (Note)

For research use only .