

产品名称: POLDIP2 Rabbit Monoclonal Antibody
产品货号: AMRe02456



产品概述 (Summary)

产品名称 (Production Name)	POLDIP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POLDIP2
别名 (Alternative Names)	p38; POLD4; PDIP38
基因 ID (Gene ID)	26073
蛋白 ID (SwissProt ID)	Q9Y2S7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 42 kDa; Observed MW: 38 kDa

研究背景 (Background)

Involved in DNA damage tolerance by regulating translesion synthesis (TLS) of templates carrying DNA damage lesions

产品名称: POLDIP2 Rabbit Monoclonal Antibody
产品货号: AMRe02456




such as 8oxoG and abasic sites.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of POLDIP2 in HeLa lysates using POLDIP2 antibody.

注意事项 (Note)

For research use only .