

产品名称: PPP2R5E Rabbit Monoclonal Antibody
产品货号: AMRe02467



产品概述 (Summary)

产品名称 (Production Name)	PPP2R5E Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP2R5E
别名 (Alternative Names)	B56E; B56epsilon
基因 ID (Gene ID)	5529
蛋白 ID (SwissProt ID)	Q16537.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 55 kDa

研究背景 (Background)

The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of

产品名称: PPP2R5E Rabbit Monoclonal Antibody
产品货号: AMRe02467




the catalytic enzyme to a particular subcellular compartment.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of PPP2R5E in C6, 3T3, mouse brain, rat brain lysates using PPP2R5E antibody.



Western blot analysis of PPP2R5E in Human placenta, Hela, 293T lysates using PPP2R5E antibody.

注意事项 (Note)

For research use only .