

产品名称: Prion Protein Rabbit Monoclonal Antibody
产品货号: AMRe02472



产品概述 (Summary)

产品名称 (Production Name)	Prion Protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRNP
别名 (Alternative Names)	CJD; GSS; PrP; PRNP; ASCR; KURU; PRIP; PrPc; CD230; AltPrP; Prion protein; PrP27-30; PrP33-35C
基因 ID (Gene ID)	5621
蛋白 ID (SwissProt ID)	P04156.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

研究背景 (Background)

产品名称: Prion Protein Rabbit Monoclonal Antibody
产品货号: AMRe02472




The PRNP gene encodes the major prion protein (PrP, CD230), a widely-expressed glycoprotein expressed at high levels in the central nervous system. While the typical cellular function of PrP is not well defined, it is a putative antioxidant and a metal-binding protein that may be involved in signal transduction.


研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Western blot analysis of Prion Protein in U251 lysates using Prion Protein antibody.



Western blot analysis of Prion protein in rat Brain and mouse Brain lysates using Prion Protein antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Prion Protein antibody. High-pressure and

产品名称: Prion Protein Rabbit Monoclonal Antibody
产品货号: AMRe02472



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .