

产品名称: Phospho-Progesterone Receptor (Ser190) Rabbit Monoclonal Antibody
产品货号: AMRe02476



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Progesterone Receptor (Ser190) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PGR
别名 (Alternative Names)	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
基因 ID (Gene ID)	5241
蛋白 ID (SwissProt ID)	P06401.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa; Observed MW: 90,118 kDa

**产品名称: Phospho-Progesterone Receptor (Ser190) Rabbit
Monoclonal Antibody**
产品货号: AMRe02476




研究背景 (Background)

Progesterone plays a central role in the reproductive events associated with the establishment and maintenance of pregnancy. Progesterone receptor, a member of the steroid receptor superfamily, mediates the physiologic effects of progesterone.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Progesterone Receptor (Phospho-Ser190) in T47D lysates using Phospho-Progesterone Receptor (Ser190) antibody.

注意事项 (Note)

For research use only .