

产品名称: Prothrombin Rabbit Monoclonal Antibody
产品货号: AMRe02483



产品概述 (Summary)

产品名称 (Production Name)	Prothrombin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	F2
别名 (Alternative Names)	Cf2; FII; Cf-2
基因 ID (Gene ID)	14061.0
蛋白 ID (SwissProt ID)	P19221.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

Thrombin, which cleaves bonds after Arg and Lys, converts fibrinogen to fibrin and activates factors V, VII, VIII, XIII, and, in

产品名称: Prothrombin Rabbit Monoclonal Antibody
产品货号: AMRe02483




complex with thrombomodulin, protein C. Functions in blood homeostasis, inflammation and wound healing .

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Prothrombin in mouse liver, rat liver lysates using Prothrombin antibody.

注意事项 (Note)

For research use only .