

产品名称: PRPF4 Rabbit Monoclonal Antibody
产品货号: AMRe02485



产品概述 (Summary)

产品名称 (Production Name)	PRPF4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRPF4
别名 (Alternative Names)	PRP4; RP70; HPRP4; Prp4p; HPRP4P; SNRNP60
基因 ID (Gene ID)	9128
蛋白 ID (SwissProt ID)	O43172.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa; Observed MW: 58 kDa

研究背景 (Background)

Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.


产品名称: PRPF4 Rabbit Monoclonal Antibody
产品货号: AMRe02485



研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of PRPF4 in K562, 3T3, HeLa lysates using PRPF4 antibody.

注意事项 (Note)

For research use only .