

产品名称: PRPF8 Rabbit Monoclonal Antibody
产品货号: AMRe02486



产品概述 (Summary)

产品名称 (Production Name)	PRPF8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRPF8
别名 (Alternative Names)	HPRP8; p220; PRP8; PRPC8; Prpf8; RP13; SNRNP220
基因 ID (Gene ID)	10594
蛋白 ID (SwissProt ID)	Q6P2Q9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 274 kDa; Observed MW: 274 kDa

研究背景 (Background)

Central component of the spliceosome, which may play a role in aligning the pre-mRNA 5'- and 3'-exons for ligation.

产品名称: PRPF8 Rabbit Monoclonal Antibody
产品货号: AMRe02486




Interacts with U5 snRNA, and with pre-mRNA 5'-splice sites in B spliceosomes and 3'-splice sites in C spliceosomes.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of PRPF8 in Jurkat, Raji lysates using PRPF8 antibody.

注意事项 (Note)

For research use only .