

产品名称: SPI1 Rabbit Monoclonal Antibody
产品货号: AMRe02497



产品概述 (Summary)

产品名称 (Production Name)	SPI1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SPI1
别名 (Alternative Names)	SPI1; Transcription factor PU.1; 31 kDa-transforming protein
基因 ID (Gene ID)	6688
蛋白 ID (SwissProt ID)	P17947.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 31 kDa; Observed MW: 42 kDa

研究背景 (Background)

Binds to the PU-box, a purine-rich DNA sequence (5'-GAGGAA-3') that can act as a lymphoid-specific enhancer. This

产品名称: SPI1 Rabbit Monoclonal Antibody
产品货号: AMRe02497




protein is a transcriptional activator that may be specifically involved in the differentiation or activation of macrophages or B-cells. Also binds RNA and may modulate pre-mRNA splicing.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of PU.1/Spi1 in Daudi, mouse liver lysates using SPI1 antibody.

注意事项 (Note)

For research use only .