

产品名称: PUMA Rabbit Monoclonal Antibody
产品货号: AMRe02498



产品概述 (Summary)

产品名称 (Production Name)	PUMA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BBC3
别名 (Alternative Names)	BBC3; PUMA; Bcl-2-binding component 3; JFY-1; p53 up-regulated modulator of apoptosis
基因 ID (Gene ID)	27113
蛋白 ID (SwissProt ID)	Q9BXH1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 21 kDa

研究背景 (Background)

产品名称: PUMA Rabbit Monoclonal Antibody
产品货号: AMRe02498




Essential mediator of p53/TP53-dependent and p53/TP53-independent apoptosis. Functions by promoting partial unfolding of BCL2L1 and dissociation of BCL2L1 from p53/TP53. Regulates ER stress-induced neuronal apoptosis.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Immunocytochemistry analysis of PUMA (green) in 293T using PUMA antibody, and DAPI (blue).




Western blot analysis of PUMA in K562 lysates using PUMA antibody.



Western blot analysis of PUMA in Jurkat, THP-1, C6, 3T3, HeLa lysates using PUMA antibody

产品名称: PUMA Rabbit Monoclonal Antibody
产品货号: AMRe02498



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PUMA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .