

产品名称: Rab45 Rabbit Monoclonal Antibody
产品货号: AMRe02506



产品概述 (Summary)

产品名称 (Production Name)	Rab45 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASEF
别名 (Alternative Names)	Ras-related protein Rab-45
基因 ID (Gene ID)	158158
蛋白 ID (SwissProt ID)	Q8IZ41.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 83 kDa; Observed MW: 83 kDa

研究背景 (Background)

Binds predominantly GDP, and also GTP.


产品名称: Rab45 Rabbit Monoclonal Antibody
产品货号: AMRe02506




研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of RAB45 in mouse heart lysates using Rab45 antibody.



Western blot analysis of RAB45 in Jurkat lysates using RAB45 antibody.



Immunocytochemistry analysis of RAB45 (green) in Hela using RAB45 antibody, and DAPI (blue)

注意事项 (Note)

产品名称: Rab45 Rabbit Monoclonal Antibody
产品货号: AMRe02506



For research use only .