

产品名称: Rab9 Rabbit Monoclonal Antibody
产品货号: AMRe02510



产品概述 (Summary)

产品名称 (Production Name)	Rab9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB9A
别名 (Alternative Names)	RAB9A; DmRab9; Ras related protein Rab 9A
基因 ID (Gene ID)	9367
蛋白 ID (SwissProt ID)	P51151.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

Involved in the transport of proteins between the endosomes and the trans Golgi network. Interacts (preferentially in its

产品名称: Rab9 Rabbit Monoclonal Antibody
产品货号: AMRe02510




GTP-bound form) with GCC2 (via its GRIP domain).


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Rab9 in 3T3, HeLa lysates using Rab9 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Rab9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .