

产品名称: RAMP1 Rabbit Monoclonal Antibody
产品货号: AMRe02519



产品概述 (Summary)

产品名称 (Production Name)	RAMP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAMP1
别名 (Alternative Names)	CRLR activity-modifying protein 1; Calcitonin-receptor-like receptor activity-modifying protein 1
基因 ID (Gene ID)	10267
蛋白 ID (SwissProt ID)	O60894.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

研究背景 (Background)

产品名称: RAMP1 Rabbit Monoclonal Antibody
产品货号: AMRe02519




Transports the calcitonin gene-related peptide type 1 receptor (CALCRL) to the plasma membrane. Acts as a receptor for calcitonin-gene-related peptide (CGRP) together with CALCRL.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of RAMP1 in Raji, COS7, C6, 3T3, HeLa lysates using RAMP1 antibody.

注意事项 (Note)

For research use only .