

产品名称: RAP1 Rabbit Monoclonal Antibody
产品货号: AMRe02521



产品概述 (Summary)

产品名称 (Production Name)	RAP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAP1B
别名 (Alternative Names)	Ras-related protein Rap-1b; K-REV; RAL1B
基因 ID (Gene ID)	5908
蛋白 ID (SwissProt ID)	P61224.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 21 kDa

研究背景 (Background)

GTP-binding protein that possesses intrinsic GTPase activity. Contributes to the polarizing activity of KRIT1 and CDH5 in the

产品名称: RAP1 Rabbit Monoclonal Antibody
产品货号: AMRe02521




establishment and maintenance of correct endothelial cell polarity and vascular lumen. Required for the localization of phosphorylated PRKCZ, PARD3 and TIAM1 to the cell junction. Plays a role in the establishment of basal endothelial barrier function.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of RAP1A/B in K562, rat Brain, C6, 3T3, HeLa lysates using RAP1 antibody.

注意事项 (Note)

For research use only .